


# AFTER 30 YEARS STILL AS EASY AS 1, 2, 3

## *1. Set Dose*



## *2. Test Impedance*



## *3. Treat*



***“THYMATRON® Why Trust Anything Less?”***




**SOMATICS, LLC**  
720 Commerce Drive, Suite 101  
Venice, FL 34292  
[www.thymatron.com](http://www.thymatron.com)

**USA toll-free 1-800-642-6761**  
(847) 234-6761  
FAX (847) 234-6763  
[sales@thymatron.com](mailto:sales@thymatron.com)

# STREAMLINED ECT: THYMATRON®

**TRY THIS**




**Bitemporal/frontal ECT**

**INSTEAD OF THIS**



**Elapsed Seizure Duration**



**Stimulus Selection**



**Titrate as % of Age\***

MECTA SPECTRUM 5000Q

| PW (msec) | F (/sec) | D (sec) | I (amp) |
|-----------|----------|---------|---------|
| 1.0       | 40       | 0.5     | 0.8     |
| 1.0       | 40       | 0.75    | 0.8     |
| 1.0       | 40       | 1.25    | 0.8     |
| 1.0       | 40       | 2.0     | 0.8     |
| 1.0       | 60       | 2.0     | 0.8     |
| 1.0       | 90       | 2.0     | 0.8     |
| 1.0       | 90       | 3.0     | 0.8     |
| 1.0       | 90       | 4.0     | 0.8     |

\*J.ECT 2013: 29:18-20